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
Pathfinder® ENDOSCOPE OVERTUBE

INDICATION FOR USE

The Pathfinder® Endoscope Overtube is intended to be used with an endoscope to facilitate intubation, change of endoscopes, and treatment in the gastrointestinal (GI) tract in adult patients (22 years of age and older).

CONTRAINDICATIONS

Use contraindicated in patients with esophageal bleeding, lesion(s), and/or laceration; esophageal strictures and/or varices; laryngeal perforation; trauma to teeth, gums, and/or pharynx; aspiration pneumonia; or any other condition that may preclude endoscopy.

DEVICE COMPATIBILITY

Refer to package labels for compatible endoscope sizes. Compatibility with other devices has not been established.

Refer to instructions for use provided with compatible devices for important use and safety information.

SPECIFICATIONS

Refer to package label(s) for device dimensions.

DESCRIPTION

Figure 1: Diagram of Pathfinder® Endoscope Overtube


Figure 1 Key

A	Handle	G	Body
B	Handle stopper	H	Tip
C	Irrigation / water line	I	Total length
D	Irrigation / water luer	J	Working length
E	Vacuum line	K	Overtube outer diameter
F	Barb	L	Handle rotator (see figure 2 on sheet 4)

POTENTIAL COMPLICATIONS

Use of the Pathfinder® Endoscope Overtube may cause trauma to the patient's GI tract.

EXPECTED CLINICAL BENEFITS

The Pathfinder® Endoscope Overtube minimizes scope looping and complications from scope looping, reduces endoscopists workload during the endoscope procedure, enhances endoscope tip stabilization and control, enables endoscopic procedure in patients with altered anatomy, enables crossing of anastomosis, and enables antegrade and retrograde enteroscopy.

WARNINGS AND PRECAUTIONS

Those who have had extensive abdominal surgeries may be poor candidates because of adhesions or altered anatomy.

The Pathfinder® Endoscope Overtube is intended for single use. **DO NOT** re-sterilize or reuse. Re-sterilization and/or reuse may compromise structural integrity of device or increase risk of contamination or infection leading to device failure and/or cross-infection and potential patient injury, illness, or death.

DO NOT use if packaging or labeling has become damaged or compromised.

DO NOT use this device for any purpose other than its intended use.

The Pathfinder® Endoscope Overtube should only be used by medical specialists with proper training in endoscopic procedures. Use proper clinical judgment regarding clinical procedures.

Operate overtube slowly to prevent damaging tissue.

Performance may be compromised if recommended lubrication is not applied to both the inner and outer surfaces of the overtube and endoscope.

The Pathfinder® Endoscope Overtube should only be used in conjunction with an endoscope. **DO NOT** insert or advance the overtube without an endoscope in place.

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DO NOT bend or pinch overtube excessively as it may cause kinking or other damage.

DO NOT use the overtube if the outer diameter of the endoscope is incompatible with the overtube. Note that the transition between the endoscope insertion tube and bending portion may have a larger outer diameter than the rest of the endoscope.

To avoid potential complications, **DO NOT** adjust overtube position if endoscope smaller than compatible endoscope size is used.

RECOMMENDED PROCEDURE

Device Preparation

1. Remove pouch from carton.
2. Use sterile technique to remove overtube from pouch.
3. Remove overtube from handle card.
4. Visually inspect device for irregularities such as sharp protrusions. Do not use device if damaged.
5. Wipe endoscope with water. Water-based lubricant (example: Medline E-Z Lubricating Jelly) may additionally be used on endoscope as desired.
6. Soak internal surface of overtube with water to activate hydrophilic coating by holding overtube in a U shape and completely filling the overtube with water for 10 seconds.
7. Empty water from overtube.
8. Add a generous amount of lubricant into handle lumen before loading overtube onto endoscope.
9. Slide endoscope through proximal end of overtube until handle of overtube contacts handle of endoscope.
10. Apply lubricant to external surface of overtube working length, overtube distal tip, and exposed portion of endoscope insertion tube.

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11. Confirm endoscope can slide smoothly through overtube. Add more water or lubricant if necessary.
12. Securely attach vacuum source to suction canister vacuum input.
13. Securely attach patient side of suction canister to overtube.
14. Turn on vacuum source (530 mmHg to 660 mmHg).
15. Refer to Figure 2. Switch overtube to rigid condition with handle rotator. Confirm overtube is rigid.
16. Switch overtube to flexible condition with handle rotator. Confirm overtube is flexible.

Endoscope and Overtube Advancement/Retraction

17. Insert endoscope into patient. Overtube should be flexible.
18. Advance/retract endoscope using standard endoscopic techniques.
19. When needed, carefully advance the flexible overtube over the endoscope without advancing the overtube beyond the endoscope tip. The flexible overtube may be advanced or retracted simultaneously with the endoscope or by alternating endoscope and overtube motion. Only advance or retract the overtube when it is in the flexible condition.
20. Switch overtube to rigid condition when needed. Do not advance, withdraw, torque, or otherwise move overtube relative to the patient when overtube is rigid. Do not move the patient when overtube is rigid.

DO NOT adjust position of the patient or the overtube relative to the patient when the overtube is rigid.

21. Advance/retract endoscope through rigidized overtube to desired location using standard endoscopic techniques. Keep rigid overtube position stationary with respect to patient. Hold overtube in place when rigid.
22. If the friction increases between the endoscope and overtube, use a syringe to insert water through the overtube irrigation luer.
23. Repeat advancement steps above as necessary to reach target location(s) and perform planned diagnostic or therapeutic procedures.

Endoscope Removal/Re-insertion

24. At any point during the procedure, it may be desirable to remove the endoscope while keeping the overtube in place within the anatomy. Overtube may be rigid or flexible during this portion of the procedure. Hold overtube in place when rigid.

DO NOT adjust position of the patient or the overtube relative to the patient when the overtube is rigid.

25. Unlock and straighten the endoscope tip.
26. Slowly withdraw the endoscope while keeping the overtube stationary. If abnormal resistance is encountered, advance the endoscope a short distance and slowly withdraw again. Do not apply excessive force to the endoscope or overtube.
27. If a new endoscope will be inserted, ensure compatibility with overtube.
28. Apply lubricant to the endoscope being inserted if necessary.
29. Hold the overtube stationary and slowly insert the endoscope through the handle of the overtube to the desired location.


Overtube Final Removal

30. At end of procedure, switch overtube to flexible condition.
31. Slowly withdraw endoscope and overtube together.
32. Disconnect overtube from vacuum source.
33. Slide overtube off endoscope.
34. Dispose of used overtube in accordance with applicable regulations.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer or its authorized representative, and the competent authority of the Member State in which the user and/or patient is established.

Symbol	Meaning	Standard	Symbol	Meaning	Standard
REF	Catalog number	ISO 15223-1:2021 Clause 5.1.6	LOT	Lot number	ISO 15223-1:2021 Clause 5.1.5
	Keep dry	ISO 15223-1:2021 Clause 5.3.4		Keep away from sunlight	ISO 15223-1:2021 Clause 5.3.2
	Use by date	ISO 15223-1:2021 Clause 5.1.4		Do not use if package has been damaged and consult instructions for use	ISO 15223-1:2021 Clause 5.2.8
	Caution	ISO 15223-1:2021 Clause 5.4.4	STERILE EO	Sterilized using ethylene oxide	ISO 15223-1:2021 Clause 5.2.3
	Do not re-sterilize	ISO 15223-1:2021 Clause 5.2.6		Manufacturer	ISO 15223-1:2021 Clause 5.1.1
	Single use only	ISO 15223-1:2021 Clause 5.4.2		Consult instructions for use	ISO 15223-1:2021 Clause 5.4.3
	Not made with natural rubber latex	USA 21 CFR 801.437	Rx Only	Caution: Federal law restricts this device to sale by or on the order of a physician	USA 21 CFR 801.109
MD	Medical device	ISO 15223-1:2021 Clause 5.7.7		Single sterile barrier system	ISO 15223-1:2021 Clause 5.2.11
CE xxx	CE Marking	Europe 2017/745		Single sterile barrier system with protective packaging outside	ISO 15223-1:2021 Clause 5.2.14
	Importer	ISO 15223-1:2021 Clause 5.1.8	OD	Compatible with endoscope insertion tube outer diameters	Not applicable
	Size	Not applicable	EU REP	Authorized representative	ISO 15223-1:2021/ Amd 1:2025 Clause 5.1.2
	Date of Manufacture	ISO 15223-1:2021 Clause 5.1.3			

Figure 2: Handle Rotator Positions

Neptune Medical, Inc.
1828 El Camino Real
Suite 508
Burlingame, CA, 94010
USA
+1 (888) 202-8401
<https://gipathfinder.com/patents/>

EU REP

Advena Ltd. Tower Business Centre, 2nd Fl.,
Tower Street, Swatar, BKR 4013 Malta
vigilance@advena.mt
Advena Services Ltd
Tower Business Centre, Tower Street
Swatar BKR 4013 Malta

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